CHANNELS[®]

MTA Reparative Cement

ENGLISH

INTRODUCTION

CHANNELS[®] MTA Reparative Cement is an endodontic cement composed of several mineral oxides. It is constituted by thin hydrophilic particles. It is indicated in cases of root canal lateral and furcation perforations, internal resorption, reverse root filling, pulp capping, pulpotomy, apexification, and apexogenesis.

MTA (Mineral Trioxide Aggregate) provides the following benefits:

- Particle size which allows complete wetting during mixing ^{1,2}.
- Excellent marginal sealing; avoids penetration of tissue fluids in the root canal ^{3,4,5}.
- Enclosing of root canal and furcation perforations through induction of periradicular cement formation ^{6,2}.
- Promotes the formation of a dentin bridge when used in pulp capping ^{3,7}.
- Unlike other cements, which demand a completely dry field, MTA is indicated when moisture control is inadequate (e.g., surgery for treatment of root perforation, reverse root filling), without loss of its properties ².

COMPOSITION

- SiO₂, K₂O, Al₂O₃, Na₂O, Fe₂O₃, SO₃, CaO, Bi₂O₃, MgO.
- Insoluble residues of CaO, KSO₄, NaSO₄ and crystalline silica.

PROPERTIES

- Setting reaction: When mixed with distilled water it forms a gel that solidifies if kept in a wet environment. The initial setting time is approximately 10 minutes and the final is 15 minutes. It is not necessary to wait for the final set to continue treatment procedures ².
- Hydrogen ion concentration (pH): After mixing, pH value is 10; in 3 hours it becomes highly alkaline (12) ^{1,8,9}.
- Radiopacity: Nearly matches that of gutta-percha. More radiopaque than dentine and bone ⁹.
- Compressive Strength: 40 MPa after 24 hours and 65 MPa after 21 days ⁹.

- ① Sites of application do not receive direct occlusal load.

DIRECTIONS FOR USE

1. Sterilize a glass slab, a metal spatula and all instruments for the insertion of CHANNELS[®] MTA Reparative Cement.
2. Mix for 30 seconds the content of 1 sachet of CHANNELS[®] MTA Reparative Cement (or 1 spoon of CHANNELS[®] MTA Reparative Cement) with 1 drop of distilled water. The mixture should be homogeneous and with a consistency similar to wet sand.
3. Place the cement on the selected site with a sterilized amalgam carrier or other appropriate instrument.
4. Condense the cement with instruments such as amalgam condensers, a number 1 spatula or

absorbent paper points moistened with distilled water.

- ① **IMPORTANT: If CHANNELS[®] MTA Reparative Cement is not used immediately after mixing, its dehydration can be prevented and the working time increased by covering the mix on the glass slab with a wet gauze. If the mixed CHANNELS[®] MTA Reparative Cement dehydrates, it must be discarded.**

INDICATIONS

1. Treatment of perforations of root canal and furcation caused iatrogenically or by caries lesion² (Pics. 1 and 2).
2. Via canal treatment of root perforation due to internal resorption² (Pic. 3).
3. Surgical treatment of root perforation due to internal resorption² (Pic. 4).
4. Periapical surgery with reverse filling⁵ (Pic. 5).
5. Pulp capping⁷.
6. Pulpotomy (removal of affected coronal pulp to preserve vitality of remaining pulp tissue² (Pic. 6).
7. Apexogenesis (induction of root development in vital teeth with an inflamed coronal pulp)².
8. Apexification (induction of formation of a mineralized barrier at the root tip of young permanent teeth with incomplete root development and a necrotic pulp² (Pic. 7).

WARNINGS AND PRECAUTIONS

- Use eye protection, mask and gloves when handling CHANNELS[®] MTA Reparative Cement. In case of eye or skin contact, wash abundantly with water.
- Only open the sachet or bottle immediately prior to use. The CHANNELS[®] MTA Reparative Cement powder is very sensitive to humidity².
- Do not use CHANNELS[®] MTA Reparative Cement to fill a root canal. Its viscosity is inadequate for this procedure and very difficult to remove in case of a reentry.
- Do not use CHANNELS[®] MTA Reparative Cement in areas of the tooth in contact with the gingival sulcus or it will be completely dissolved.
- Only use CHANNELS[®] MTA Reparative Cement after remission of acute signs and symptoms of the endodontic disease. The acidic pH of endodontically compromised sites (lesions) prevents its setting reaction.
- Apply CHANNELS[®] MTA Reparative Cement carefully. Similar to other endodontic cements, it is resorbed if extruded. However, the excess of any cement may impair the healing process.
- ① **IMPORTANT: The information provided in this manual is based on laboratory and clinical studies. The successful use of CHANNELS[®] MTA Reparative Cement depends on a correct diagnosis, the operative technique, the condition of the treated tooth and the general health of the patient. This product must be used according to this manual.**

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34:80-82.

5. Torabinejad M, Hong C, Lee SJ, Monsef M, Pitt Ford TR. Investigation of mineral trioxide aggregate for root-end filling in dogs. J Endod. 1995 Dec; 21(12):603-608.
6. Torabinejad, M., Pitt Ford, T., McKentry, D., Abedi, H., Miller, D., Kariyavasam, S. Histologic assessment of mineral trioxide aggregate as a root-end filling in monkeys; J Endod; 1997; 23:225-228.
7. Faraco JR IM, Holland R. Response of the pulp of dogs to capping with mineral trioxide aggregate or a calcium hydroxide cement. Dent Traumatol. 2001 Aug; 17(4):163-6.
8. Parirokh, M. and Torabinejad, M. Mineral Trioxide Aggregate: A Comprehensive Literature Review - Part III: Clinical Applications, Drawbacks, and Mechanism of Action; J Endod; 2010; 36:400-413
9. Torabinejad M, Hong C, McDonald F, Pitt Ford TR. Physical and chemical properties of a new root endfilling material. J Endod. 1995 Jul; 21(7):349-53.
10. Rodrigo Ricci Vivian, et al. Evaluation of the physical and chemical properties of two commercial and three experimental root-end filling materials; Oral Surgery, Oral Medicine, Oral Pathology, Oral Radiology and Endodontology, 2010; 110:433-441.

RX ONLY

Manufactured for:
Insight Endo
41 Weaver Road
Denver, PA 17517 USA

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ATTENTION: This product must be used according to the instructions described in this manual. The manufacturer is not responsible for failure or damage caused by incorrect handling or use.

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